DFTA – MONTHLY ONSITE AED INSPECTION FORM

Date:				
AED Unit Clean, no dirt or contamination; no damage present				
 Supplies Available: a. Pads A: Installed in unit within expiration date undamaged b. Pads B: Stored in lid of carrying case within expiration date undamaged c. Ancillary Supplies – Fast Response Kit -Red pouch contains: <i>CPR mask</i> <i>Gloves</i> <i>Scissors</i> <i>Razor</i> <i>Towel</i> d. Spare battery stored in lid of carrying case - unopened battery within install before date (optional for new AED) 				
Status Indicator green flashing light.				
Inspected By:				
Remarks, Problems, Corrective Actions:			 	

Prepared by Emergency Skills, Inc. Daily, check AED status indicator. Monthly, perform thorough check, Date and initial upon completion. Please make copies of this form as needed. Store completed forms in a safe place.

To order supplies, call Emergency Skills, Inc. at 212-564-6833